

SSCOR

SSCOR VX-2

Complies with SAE J3043 Standard!



SSCOR VX-2

The SSCOR VX-2® is a portable, powerful battery powered suction unit that can be mounted in the action area of an ambulance and used as onboard suction. The optional charging/retention bracket connects the VX-2 to vehicle power at all times. With its variable regulator and gauge, outward facing controls, and drawing its power from the vehicle, the VX-2 is setup for on-board

suction. If portable suction is required at the scene, the VX-2 can easily be removed from the retention bracket and used as portable suction, drawing power from its battery.

The VX-2 comes equipped with a battery care system which insures the battery is charged to an optimal level, prevents deep discharge and indicates battery condition.

Complete information at www.sscor.com.



Use the SSCOR Onboard Canister Clip #230-00003 (sold separately) and a Bemis 1200cc canister as shown above for SAE J3043 compliance.



SAE J3043 and EN1789 compliant Charging/Retention Bracket. Holds to a force of 25 G and conducts vehicle power to run and charge the unit.

OPTIONS!

Different solutions for different needs.



	VX-2 Model #2310BV	VX-2 Model #2310V
Size	17" L x 9" H x 5.25" W 43.2cm L x 22.9cm H x 13.3cm W	
Weight	8.6 lbs. or 3.9 kg	
Vacuum Pump	12V DC oil-less diaphragm. Clinical Airflow >=30LPM	
Regulator	Variable Regulator and Gauge	
Canister	Bemis canister. For canister options, see http://www.sscor.com/S-SCORT_VX2_canister_battery.html	
Run Time (full capacity)	45 minutes (+/- 10%)	
Battery Condition Indication	3 LED, 5 stage, battery condition indicator system	
Charging	Charge constantly from vehicle via the charging retention bracket or an optional AC/DC converter Part # 80521-100 for land based AC/DC charging	Charge constantly from vehicle power using DC power cord (included) or an AC/DC converter Part # 80521-100 for land based AC/DC charging
Charging/Retention Bracket	Part # 8310 Included SAE J3043 Compliant, EN1789 Compliant	No
DC Power Cord	Optional	Part # 80665
230 Setup International Customers	Optional AC/DC charger Part # 80529	DC Power Cord Part # 80965; AC/DC charger Part # 80529
Items Included	Suction Unit, Charging/Retention Bracket, 1200cc Bemis® Canister, patient tubing and HI-D® suction tip	Suction Unit, DC Power Cord, AC/DC Converter, 1200cc Bemis Canister, patient tubing and HI-D suction tip
Onboard Canister Clip	Optional Part #230-00003	N/A



The NEW Onboard Canister Clip

Use the SSSCOR Onboard Canister Clip with the VX-2 Model 2310BV and the Model 23002 Canister Holder to comply with SAE J3043. The Onboard Canister Clip will hold the Bemis 1200cc canister in place on the VX-2 Model 2310BV and 23002 plastic canister holder to a force of 25G.
Part # 230-00003