



INTRODUCING THE S-SCORT® S3™ OROPHARYNGEAL EVACUATION TOOL


The S-SCORT S3 is a rigid yet flexible large bore suction tip with a concave "scooping" surface which allows for removal of solids and other material that will not travel through a large bore suction tube.



By placing the slider over the thumb port, the S3 can be briefly placed in the deep oropharynx to suction, freeing both hands to manage intubation.



The integral stylet allows the S3 to flex and hold its' shape!!




US PATENT PENDING

Order Number: 44305

For more Information, please contact:



LIT-44305 Rev. A



Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: (+31) 70 345 8570
Fax: (+31) 70 346 7299



11064 Randall Street • Sun Valley, CA • 91352 USA
818-504-4054 • 800-434-5211
marketing@sscor.com • www.sscor.com