

Power In The Palm Of Your Hand

Now with a 400cc Canister!



The Quickdraw

The **SSCOR** Quickdraw easily fits into an airway or first-in bag so you can have quality medical suction available at the patient's side. The **SSCOR** Quickdraw provides suction as efficient as a large aspirator but weighs only 2.45 pounds (1.10kg).



NEVER NEEDS CHARGING!

Powered by alkaline battery technology!

Alkaline battery technology is ideal for small portable devices carried in kits. Alkaline batteries have a long shelf life, require no charging, pack a lot of power for long run times and are widely available and inexpensive.

SSCOR offers an alkaline battery pack which will power the Quickdraw for 3 hours of continuous use. And now **SSCOR** also offers the all new 'AAA' battery holder which accepts standard off the shelf 'AAA' alkaline batteries.



The Quickdraw is the ideal aspirator to carry in your 'first in' kit and have on hand every time you are with the patient.

Quickdraw! Have it ready to use at your side on every call!



Items Included



400cc Canister with Hose Barb

The Quickdraw canister with hose barb will accept standard suction tubing. For use with common suction tips including tracheal catheters.

#2488 Each

#2488-20 Case of 20



SSCOR SDC Suction Catheter™

(Formerly the SSCOR DuCanto Catheter.) A large bore suction catheter (.26" I.D.) designed for use during routine and emergency airway management.

#200-00002 Each

#200-00002C Case of 50



Battery Pack

Run time for the SSCOR alkaline battery pack is three hours +/- 10%.

#80613-100 Each



Connecting Tubing 3 Feet

3' suction connecting tubing.

#43203 Each

#43203C Case of 50



Battery Holder

Accepts 10 'AAA' batteries. Run times with a SSCOR battery holder will vary from 60 to 100 minutes depending on the brand and freshness of batteries used. (Batteries not included.)

#80615-100 Each



Optional Products



SSCOR HI-D "Big Stick" Suction Tip

A large bore suction tip (.25" I.D.) with thumb port, allowing the patients airway to be cleared more effectively.

#44241 Each

#44241C Case of 50

Specifications

Size:

10.5" L x 4.5" H x 4.35" W (27cm L x 11cm H x 11cm)

Weight:

2.45 lbs or 1.10 kg

Performance:

Maximum negative pressure >500mmHg; Low negative pressure setting approximately 80 - 100mmHg

Collection Canister:

#2488 400 cc/ml capacity; 99.98% bacterial effective filter

Run Time:

#80613-100 battery pack: 3 hours +/- 10%

#80615-100 battery holder: 60 - 100 minutes with 10 new 'AAA' batteries (batteries not included)

Visit shop.sscor.com to purchase the Quickdraw, replacement parts, disposables and accessories. Complete product information can be found at www.sscor.com.

SSCOR, Inc.

+1 818 504-4054

2125 N. Madera Rd. Unit C, Simi Valley, CA 93065 USA

e-mail: sales@sscor.com www.sscor.com shop.sscor.com

© 2022, Quickdraw is a registered trademark of SSCOR, Inc., Patents 7,063,688 B2 - D564,654 S - 7,938,794 B2 LIT-2403-1